
Fluoromount-G[®] Anti-Fade Mounting Medium




Quality Antibodies for Quality Research



Mounting Medium	Cat. No.	Qty
Fluoromount-G [®] Anti-Fade	0100-35	20 mL

Fluoromount-G[®] Anti-Fade is an aqueous-based compound for mounting slides used in microscopy applications in which the preservation of fluorescence intensity is critical. This mounting medium exhibits minimal background in the blue spectrum and provides a semi-permanent seal for storage of slide preparations.



Corporate Offices

160 Oxmoor Blvd Birmingham, AL 35209 USA


Mailing Address

PO Box 26221 Birmingham, AL 35260 USA


Tel 205.945.1774 Fax 205.945.8768 Email info@SouthernBiotech.com Technical Support technical@SouthernBiotech.com Website www.SouthernBiotech.com

SouthernBiotech is ISO 9001:2015 Certified





Frozen newborn mouse cartilage section was blocked with Normal Rabbit Serum (SB Cat. No. 0040-01) and stained with Mouse Anti-Collagen Type II-UNLB followed by Rabbit F(ab')2 Anti-Mouse IgG(H+L), Human ads-BIOT (SB Cat. No. 6125-08), Streptavidin-FITC (SB Cat. No. 7100-02), DAPI, and mounted with Fluoromount-G[®] Anti-Fade (SB Cat. No. 0100-35).



Paraffin embedded mouse small intestine tissue section was blocked with Normal Rabbit Serum (SB Cat. No. 0040-01) and stained with Goat Anti-Mouse IgA-UNLB (SB Cat. No. 1040-01) followed by Rabbit Anti-Goat IgG(H+L), Human SP ads-BIOT (SB Cat. No. 6164-08), Streptavidin-CY3.5 (SB Cat. No. 7100-24), DAPI, and mounted with Fluoromount-G[®] Anti-Fade (SB Cat. No. 0100-35).