

Swine IgG

<u>Cat. No.</u>	<u>Form</u>	<u>Quantity</u>
0137-01	Purified (UNLB) IgG	10.0 mg

DESCRIPTION

Source Normal swine serum
Ig Isotype Swine IgG
Specificity N/A

RESEARCH APPLICATIONS

- Enzyme-Linked-Immunosorbent-Assay (ELISA)
- Flow cytometry

CHARACTERIZATION

To insure lot-to-lot consistency, each batch of product is tested by ELISA and/or flow cytometry for conformance to characteristics of a standard reference reagent.

WORKING DILUTIONS

Since applications vary, you should determine the optimum working dilution of the product that is appropriate for your specific need.

HANDLING AND STORAGE

- The purified (UNLB) IgG is supplied as 10.0 mg of purified immunoglobulin in 2.0 mL of 100 mM borate buffered saline, pH 8.2. *No preservatives or amine-containing buffer salts added.* Store at 2-8°C.

For Research Use Only. Not for Diagnostic or Therapeutic Use.