




Rat Anti-Mouse CD4

Cat. No.	Format	Size
1540-01	Purified (UNLB)	0.5 mg
1540-02	Fluorescein (FITC)	0.5 mg
1540-02S	Fluorescein (FITC)	0.1 mg
1540-08	Biotin (BIOT)	0.5 mg
1540-09	R-phycoerythrin (PE)	0.1 mg
1540-09L	R-phycoerythrin (PE)	0.2 mg
1540-10	R-phycoerythrin-Texas Red® (PE/TXRD)	0.1 mg
1540-11	Allophycocyanin (APC)	0.1 mg
1540-13	Spectral Red® (SPRD)	0.1 mg
1540-14	Low Endotoxin, Azide-Free (LE/AF)	0.5 mg
1540-15	Cyanine 5 (CY5)	0.1 mg
1540-16	R-phycoerythrin-Cyanine 5.5 (PE/CY5.5)	0.1 mg
1540-17	R-phycoerythrin-Cyanine 7 (PE/CY7)	0.1 mg
1540-18	Allophycocyanin-Cyanine 5.5 (APC/CY5.5)	0.1 mg
1540-19	Allophycocyanin-Cyanine 7 (APC/CY7)	0.1 mg
1540-26	Pacific Blue™ (PACBLU)	0.1 mg
1540-27	Alexa Fluor® 700 (AF700)	0.1 mg
1540-30	Alexa Fluor® 488 (AF488)	0.1 mg
1540-31	Alexa Fluor® 647 (AF647)	0.1 mg



BALB/c mouse splenocytes were stained with Rat Anti-Mouse CD4-AF647 (SB Cat. 1540-31) and Rat Anti-Mouse CD8α-FITC (SB Cat. No. 1550-02).

Overview

Clone	GK1.5
Isotype	Rat (Lewis) IgG _{2b} K
Immunogen	Mouse CTL clone V4
Specificity	Mouse CD4; Mr 55 kDa
Alternate Name(s)	L3T4, Ly-4

Description

CD4 is a type I transmembrane glycoprotein expressed on most thymocytes and the helper/inducer subpopulation of mature T lymphocytes of all mouse strains tested. CD4 has also been detected at low density on pluripotent hematopoietic stem cells, bone marrow myeloid precursors, and intrathymic lymphoid precursors. The CD4 antigen functions as a co-receptor, interacting with class II major histocompatibility complex (MHC) molecules in the recognition of foreign antigens by T cells. *In vivo* and *in vitro* treatment with the GK1.5 monoclonal antibody effectively depletes CD4⁺ cells. GK1.5 also blocks helper T cell responses to MHC class II antigens, including cytolysis, proliferation, allogeneic B cell help, and release of lymphokines.

Applications

FC – Quality tested ^{1,5,7,11-15}
 IHC-FS – Reported in literature ²⁻⁴
 ICC – Reported in literature ⁵
 IP – Reported in literature ¹
 Depletion – Reported in literature ^{1,7,8}
 Sep – Reported in literature ^{6,9,10}
 Block – Reported in literature ^{1,8}

Working Dilutions

Flow Cytometry	FITC, BIOT, AF488, and PACBLU conjugates	≤ 1 µg/10 ⁶ cells
	PE, APC, PE/TXRD, SPRD, CY5, PE/CY5.5, PE/CY7, APC/CY5.5, APC/CY7, AF647, and AF700 conjugates	≤ 0.2 µg/10 ⁶ cells
	For flow cytometry, the suggested use of these reagents is in a final volume of 100 µL	

Other Applications Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.

For Research Use Only. Not for Diagnostic or Therapeutic Use.

Handling and Storage

- The purified (UNLB) antibody is supplied as 0.5 mg of purified immunoglobulin in 1.0 mL of borate buffered saline, pH 8.2. *No preservatives or amine-containing buffer salts added.* Store at 2-8°C.
- The fluorescein (FITC) conjugate is supplied as 0.5 mg or 0.1 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The biotin (BIOT) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The R-phycoerythrin (PE) conjugate is supplied as 0.1 mg in 1.0 mL or 0.2 mg in 2.0 mL of PBS/NaN₃ and a stabilizing agent. Store at 2-8°C. **Do not freeze!**
- The allophycocyanin (APC) conjugate is supplied as 0.1 mg in 1.0 mL of PBS/NaN₃ and a stabilizing agent. Store at 2-8°C. **Do not freeze!**
- The Spectral Red® (SPRD), R-phycoerythrin-Texas Red® (PE/TXRD), R-phycoerythrin-Cyanine 7 (PE/CY7), R-phycoerythrin-Cyanine 5.5 (PE/CY5.5), allophycocyanin-Cyanine 5.5 (APC/CY5.5), and allophycocyanin-Cyanine 7 (APC/CY7) conjugates are supplied as 0.1 mg in 1.0 mL of PBS/NaN₃ and a stabilizing agent. Store at 2-8°C. **Do not freeze!**
- The low endotoxin, azide-free (LE/AF) antibody is supplied as 0.5 mg purified immunoglobulin in 1.0 mL of PBS. Contains no preservative; handle under aseptic conditions. Store at 2-8°C or aliquot into smaller volumes and store at -20°C. Avoid multiple freeze / thaw cycles.
- The Cyanine 5 (CY5) conjugate is supplied as 0.1 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The Alexa Fluor® 488 (AF488), Alexa Fluor® 647 (AF647), Alexa Fluor® 700 (AF700), and Pacific Blue™ (PACBLU) conjugates are supplied as 0.1 mg in 0.2 mL of PBS/NaN₃. Store at 2-8°C.
- Protect fluorochrome-conjugated forms from light. Reagents are stable for the period shown on the label if stored as directed.

Warning

Some reagents contain sodium azide. Please refer to product specific SDS.

References

1. Dialynas DP, Quan ZS, Wall KA, Pierres A, Quintáns J, Loken MR, et al. Characterization of the murine T cell surface molecule, designated L3T4, identified by monoclonal antibody GK1.5: similarity of L3T4 to the human Leu-3/T4 molecule. *J Immunol.* 1983;131:2445-51. (Immunogen, FC, IP, Depletion, Block)
2. Delong P, Tanaka T, Kruklytis R, Henry AC, Kapoor V, Kaiser LR, et al. Use of cyclooxygenase-2 inhibition to enhance the efficacy of immunotherapy. *Cancer Res.* 2003;63:7845-52. (IHC-FS)
3. Xia RH, Yosef N, Ubogu EE. Clinical, electrophysiological and pathologic correlations in a severe murine experimental autoimmune neuritis model of Guillain-Barré syndrome. *J Neuroimmunol.* 2010;219:54-63. (IHC-FS)
4. Kosuge H, Ishihara T, Haraguchi G, Maejima Y, Okada H, Saiki H, et al. Treatment with telmisartan attenuates graft arteriosclerosis in murine cardiac allografts. *J Heart Lung Transplant.* 2010;29:562-7. (IHC-FS)
5. Choi J, Oh S, Lee D, Oh HJ, Park JY, Lee SB, et al. Mst1-FoxO signaling protects naïve T lymphocytes from cellular oxidative stress in mice. *PLoS One.* 2009;4(11):e8011. (ICC, FC)
6. Alari-Pahissa E, Vega-Ramos J, Zhang J, Castaño AR, Turley SJ, Villadangos JA, et al. Differential effect of CD69 targeting on bystander and antigen-specific T cell proliferation. *J Leukoc Biol.* 2012;92:145-58. (Sep)
7. Grodeland G, Mjaaland S, Roux KH, Fredriksen AB, Bogen B. DNA vaccine that targets hemagglutinin to MHC class II molecules rapidly induces antibody-mediated protection against influenza. *J Immunol.* 2013;191:3221-31. (Depletion, FC)
8. Ye X, Robinson MB, Pabin C, Batschaw ML, Wilson JM. Transient depletion of CD4 lymphocyte improves efficacy of repeated administration of recombinant adenovirus in the ornithine transcarbamylase deficient sparse fur mouse. *Gene Ther.* 2000;7:1761-7. (Depletion, Block)
9. Travers H, Anderson G, Gentle D, Jenkinson E, Girdlestone J. Protocols for high efficiency, stage-specific retroviral transduction of murine fetal thymocytes and thymic epithelial cells. *J Immunol Methods.* 2001;253:209-22. (Sep)
10. Vega-Ramos J, Alari-Pahissa E, del Valle J, Carrasco-Marín E, Esplugues E, Borràs M, et al. CD69 limits early inflammatory diseases associated with immune response to Listeria monocytogenes infection. *Immunol Cell Biol.* 2010;88:707-15. (Sep)
11. Avagyan S, Glouchkova L, Choi J, Snoeck H. A quantitative trait locus on chromosome 4 affects cycling of hematopoietic stem and progenitor cells through regulation of TGF-β2 responsiveness. *J Immunol.* 2008;181:5904-11. (FC)
12. Barnes MJ, Krebs P, Harris N, Eidenschenk C, Gonzalez-Quintal R, Arnold CN, et al. Commitment to the regulatory T cell lineage requires CARMA1 in the thymus but not in the periphery. *PLoS Biol.* 2009;7(3):e1000051. (FC)
13. Safari D, Dekker HA, Rijkers G, Snippe H. Codelivery of adjuvants at the primary immunization site is essential for evoking a robust immune response to neoglycoconjugates. *Vaccine.* 2011;29:849-54. (FC)
14. Berge T, Grønningssæter IH, Lorvik KB, Abrahamson G, Granum S, Sundvold-Gjerstad V, et al. SH2D2A modulates T cell mediated protection to a B cell derived tumor in transgenic mice. *PLoS One.* 2012;7(10):e48239. (FC)
15. Jacobsen J, Haabeth OW, Tveita AA, Schjetne KW, Munthe LA, Bogen B. Naive idiotope-specific B and T cells collaborate efficiently in the absence of dendritic cells. *J Immunol.* 2014;192:4174-83. (FC)

Texas Red® is a registered trademark of Molecular Probes, Inc.

Spectral Red® is a registered trademark of Southern Biotechnology Associates, Inc.

Spectral Red® is a PE/CY5 tandem conjugate.

Cy™ is a trademark of Cytiva or one of its subsidiaries.

Alexa Fluor® 488, 647, 700 and Pacific Blue™ are provided under an agreement between Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corporation), and Southern Biotechnology Associates, Inc., and the manufacture, use, sale or import of this product may be subject to one or more U.S. patents, pending applications, and corresponding non-U.S. equivalents, owned by Molecular Probes, Inc. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For information on purchasing a license to this product for any other use, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, USA, Tel: (541) 465-8300. Fax: (541) 335-0504.