




## Mouse Anti-Mouse I-E<sup>k</sup>

Cat. No.	Format	Size
1897-01	Purified (UNLB)	0.5 mg
1897-02	Fluorescein (FITC)	0.5 mg
1897-08	Biotin (BIOT)	0.5 mg
1897-09	R-phycoerythrin (PE)	0.1 mg
1897-14	Low Endotoxin, Azide-Free (LE/AF)	0.5 mg



AKR mouse splenocytes were stained with Mouse Anti-Mouse I-E<sup>k</sup>-PE (SB Cat. No. 1897-09) and Rat Anti-Mouse CD3ε-FITC (SB Cat. No. 1535-02).

### Overview

<b>Clone</b>	14-4-4S
<b>Isotype</b>	Mouse (C3H.SW) IgG <sub>2aκ</sub>
<b>Immunogen</b>	C3H mouse skin graft and splenocytes
<b>Specificity</b>	Mouse I-E <sup>k</sup> /Rat RT1D
<b>Alternate Name(s)</b>	MHC Class II

### Description

The monoclonal antibody 14-4-4S reacts with the I-E<sup>k</sup> class II alloantigen on cells from mice of the H-2<sup>d</sup>, H-2<sup>p</sup>, and H-2<sup>r</sup> haplotypes. Cells from mice of the H-2<sup>b</sup>, H-2<sup>f</sup>, H-2<sup>q</sup>, and H-2<sup>s</sup> haplotypes do not express I-E antigen. The antibody has been reported to cross-react with the rat MHC class II alloantigen RT1D.

### Applications

FC – Quality tested <sup>8,9</sup>  
 IHC-FS – Reported in literature <sup>7</sup>  
 ICC – Reported in literature <sup>6</sup>  
 IP – Reported in literature <sup>2,3</sup>  
 CMCD – Reported in literature <sup>1</sup>  
 Activ – Reported in literature <sup>4</sup>  
 Block – Reported in literature <sup>5</sup>

### Working Dilutions

<b>Flow Cytometry</b>	FITC and BIOT conjugates	≤ 1 μg/10 <sup>6</sup> cells
	PE conjugate	≤ 0.2 μg/10 <sup>6</sup> cells
For flow cytometry, the suggested use of these reagents is in a final volume of 100 μL		

**Other Applications** Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.

**For Research Use Only. Not for Diagnostic or Therapeutic Use.**

## Handling and Storage

---

- The purified (UNLB) antibody is supplied as 0.5 mg of purified immunoglobulin in 1.0 mL of borate buffered saline, pH 8.2. *No preservatives or amine-containing buffer salts added.* Store at 2-8°C.
- The fluorescein (FITC) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN<sub>3</sub>. Store at 2-8°C.
- The biotin (BIOT) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN<sub>3</sub>. Store at 2-8°C.
- The R-phycoerythrin (PE) conjugate is supplied as 0.1 mg in 1.0 mL of PBS/NaN<sub>3</sub> and a stabilizing agent. Store at 2-8°C. **Do not freeze!**
- The low endotoxin, azide-free (LE/AF) antibody is supplied as 0.5 mg purified immunoglobulin in 1.0 mL of PBS. Contains no preservative; handle under aseptic conditions. Store at 2-8°C or aliquot into smaller volumes and store at -20°C. Avoid multiple freeze / thaw cycles.
- Protect fluorochrome-conjugated forms from light. Reagents are stable for the period shown on the label if stored as directed.

## Warning

---

Some reagents contain sodium azide. Please refer to product specific SDS.

## References

---

1. Ozato K, Mayer N, Sachs DH. Hybridoma cell lines secreting monoclonal antibodies to mouse H-2 and Ia antigens. *J Immunol.* 1980;124:533-40. (Immunogen, CMCD)
2. Neiß U, Reske K. Non-coordinate synthesis of MHC class II proteins and invariant chains by epidermal Langerhans cells derived from short-term in vitro culture. *Int Immunol.* 1994;6:61-71. (IP, Rat RT1D Reactivity)
3. Bhattacharya A, Dorf ME, Springer TA. A shared alloantigenic determinant on Ia antigens encoded by the I-A and I-E subregions: evidence for I region gene duplication. *J Immunol.* 1981;127:2488-95. (IP)
4. Bishop GA, Haughton G. Induced differentiation of a transformed clone of Ly-1<sup>+</sup> B cells by clonal T cells and antigen. *Proc Natl Acad Sci USA.* 1986;83:7410-4. (Activ)
5. Roy M, Aruffo A, Ledbetter J, Linsley P, Kehry M, Noelle R. Studies on the interdependence of gp39 and B7 expression and function during antigen-specific immune responses. *Eur J Immunol.* 1995;25:596-603. (Block)
6. Lenz A, Heuffer C, Rammensee H, Glassl H, Koch F, Romani N, et al. Murine epidermal Langerhans cells express significant amounts of class I major histocompatibility complex antigens. *Proc Natl Acad Sci USA.* 1989;86:7527-31. (ICC)
7. Farr A, DeRoos PC, Eastman S, Rudensky AY. Differential expression of CLIP:MHC class II and conventional endogenous peptide:MHC class II complexes by thymic epithelial cells and peripheral antigen-presenting cells. *Eur J Immunol.* 1996;26:3185-93. (IHC-FS)
8. Rinderknecht CH, Roh S, Pashine A, Belmares MP, Patil NS, Lu N, et al. DM influences the abundance of major histocompatibility complex class II alleles with low affinity for class II-associated invariant chain peptides via multiple mechanisms. *Immunology.* 2010;131:18-32. (FC)
9. Rinderknecht CH, Lu N, Crespo O, Truong P, Hou T, Wang N, et al. I-A<sup>g7</sup> is subject to post-translational chaperoning by CLIP. *Int Immunol.* 2010;22:705-16. (FC)

TB1897  
08-Oct-21