




Rat Anti-Mouse CD107b

Cat. No.	Format	Size
1921-01	Purified (UNLB)	0.5 mg
1921-02	Fluorescein (FITC)	0.5 mg
1921-11	Allophycocyanin (APC)	0.1 mg



Mouse M ϕ cells were stained with Rat Anti-Mouse CD107b-APC (SB Cat. No. 1921-11).

Overview

Clone	ABL-93
Isotype	Rat (Sprague Dawley) IgG _{2a} κ
Immunogen	Glycoproteins purified from BALB/c mouse embryo 3T3 cell line
Specificity	Mouse CD107b; Mr 100-110 kDa
Alternate Name(s)	LAMP-2, lysosomal-associated membrane protein 2

Description

CD107b, also known as lysosomal-associated membrane protein 2 (LAMP-2), is a heavily glycosylated type I transmembrane protein that constitutes one of the two major sialoglycoproteins on lysosomal membranes that can be used to distinguish lysosomes from other organelles. CD107b is similar in structure to Mac-1 and is expressed primarily on epithelial cells and macrophages in a variety of organs in normal and beige mutant mice.

Applications

FC – Quality tested ⁷
 IHC-FS – Reported in literature ⁴
 ICC – Reported in literature ^{1,2,5,6}
 EM – Reported in literature ¹
 IP – Reported in literature ^{1,2}
 WB – Reported in literature ^{2,3}

Working Dilutions

Flow Cytometry	Purified (UNLB) antibody	$\leq 1 \mu\text{g}/10^6$ cells
	FITC and APC conjugates	$\leq 1 \mu\text{g}/10^6$ cells
For flow cytometry, the suggested use of these reagents is in a final volume of 100 μL .		

Other Applications Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.

For Research Use Only. Not for Diagnostic or Therapeutic Use.

Handling and Storage

- The purified (UNLB) antibody is supplied as 0.5 mg of purified immunoglobulin in 1.0 mL of borate buffered saline, pH 8.2. *No preservatives or amine-containing buffer salts added.* Store at 2-8°C.
- The fluorescein (FITC) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The allophycocyanin (APC) conjugate is supplied as 0.1 mg in 1.0 mL of PBS/NaN₃ and a stabilizing agent. Store at 2-8°C. **Do not freeze!**
- Protect fluorochrome-conjugated forms from light. Reagents are stable for the period shown on the label if stored as directed.

Warning

Some reagents contain sodium azide. Please refer to product specific SDS.

References

1. Chen JW, Murphy TL, Willingham MC, Pastan I, August JT. Identification of two lysosomal membrane glycoproteins. *J Cell Biol.* 1985;101:85-95. (Immunogen, ICC, EM, IP)
2. Andrejewski N, Punnonen E, Guhde G, Tanaka Y, Lüllmann-Rauch R, Hartmann D, et al. Normal lysosomal morphology and function in LAMP-1-deficient mice. *J Biol Chem.* 1999;274:12692-701. (IP, WB, ICC)
3. Via LE, Deretic D, Ulmer RJ, Hibler NS, Huber LA, Deretic V. Arrest of mycobacterial phagosome maturation is caused by a block in vesicle fusion between stages controlled by rab5 and rab7. *J Biol Chem.* 1997;272:13326-31. (WB)
4. Zhao W, Dumanis SB, Tamboli IY, Rodriguez GA, LaDu MJ, Moussa CE, et al. Human APOE genotype affects intraneuronal A β ₁₋₄₂ accumulation in a lentiviral gene transfer model. *Hum Mol Genet.* 2014;23:1365-75. (IHC-FS)
5. Martina JA, Moriyama K, Bonifacino JS. BLOC-3, a protein complex containing the Hermansky-Pudlak syndrome gene products HPS1 and HPS4. *J Biol Chem.* 2003;278:29376-84. (ICC)
6. Valk E, Leung R, Kang H, Kaneko K, Rudd CE, Schneider H. T cell receptor-interacting molecule acts as a chaperone to modulate surface expression of the CTLA-4 coreceptor. *Immunity.* 2006;25:807-21. (ICC)
7. Ramachandra L, Sramkoski RM, Canaday DH, Boom WH, Harding CV. Flow analysis of MHC molecules and other membrane proteins in isolated phagosomes. *J Immunol Methods.* 1998;213:53-71. (FC)