
SouthernBiotech



Rat IgG_{2b} Isotype Control

Cat. No.	Format	Size
0127-01	Purified (UNLB)	0.5 mg



ELISA plate was coated with serially diluted Rat IgG2b-UNLB (SB Cat. No. 0127-01). Immunoglobulin was detected with Mouse Anti-Rat IgG1-HRP (SB Cat. No. 3060-05), Mouse Anti-Rat IgG2a-HRP (SB Cat. No. 3065-05), Mouse Anti-Rat IgG2b-HRP (SB Cat. No. 3070-05), Mouse Anti-Rat IgG2c-HRP (SB Cat. No. 3075-05), and Goat Anti-Rat IgM-HRP (SB Cat. No. 3030-05) and quantified.

Overview

Applications

ELISA - Quality tested

Working Dilutions

ELISA Purified (UNLB) antibody ≤ 1 μg/mL

Other Applications Since applications vary, you should determine the optimum working dilution for the product that is

appropriate for your specific need.

Handling and Storage

- The purified (UNLB) antibody is supplied as 0.5 mg of purified immunoglobulin in 1.0 mL of borate buffered saline, pH 8.2. No preservatives or amine-containing buffer salts added. Store at 2-8°C.
- · Reagent is stable for the period shown on the label if stored as directed.

For Research Use Only. Not for Diagnostic or Therapeutic Use.

TB0127 07<u>-Feb-18</u>