



Swine IgG Isotype Control

Cat. No.	Format	Size
0137-01	Purified (UNLB)	10.0 mg

Description

Isotype	Swine IgG
Source	Normal swine serum

Applications

Quality tested applications include –

FC
ELISA

Working Dilutions

ELISA	Purified (UNLB) antibody	$\leq 2 \mu\text{g/mL}$
--------------	--------------------------	-------------------------

Other Applications	Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.
---------------------------	--

Handling and Storage

- The purified (UNLB) IgG is supplied as 10.0 mg of purified immunoglobulin in 2.0 mL of borate buffered saline, pH 8.2. *No preservatives or amine-containing buffer salts added.* Store at 2-8°C.
- Reagents are stable for the period shown on the label if stored as directed.

For Research Use Only. Not for Diagnostic or Therapeutic Use.

TB0137
20-Mar-15

Corporate Offices: 160 Oxmoor Blvd • Birmingham, AL 35209 • USA **Mailing Address:** P.O. Box 26221 • Birmingham, AL 35260 • USA

Tel: 205.945.1774 • U.S. and Canada: 800.722.2255 • **Fax:** 205.945.8768

Email: info@southernbiotech.com • **Website:** www.southernbiotech.com