




Rat Anti-Mouse CD11a

Cat. No.	Format	Size
1555-01	Purified (UNLB)	0.5 mg
1555-02	Fluorescein (FITC)	0.5 mg
1555-08	Biotin (BIOT)	0.5 mg
1555-09	R-phycoerythrin (PE)	0.1 mg
1555-09L	R-phycoerythrin (PE)	0.2 mg
1555-14	Low Endotoxin, Azide-Free (LE/AF)	0.5 mg
1555-17	R-phycoerythrin-Cyanine 7 (PE/CY7)	0.1 mg



BALB/c mouse splenocytes were stained with Rat Anti-Mouse CD11a-PE/CY7 (SB Cat. No. 1555-17).

Overview

Clone	I21/7
Isotype	Rat (Lewis) IgG _{2a} κ
Immunogen	Whole BW5147 cells
Specificity	Mouse/Human CD11a; Mr 180 kDa
Alternate Name(s)	LFA-1, LFA-1α, integrin α _L

Description

The monoclonal antibody I21/7 reacts with the α subunit of mouse and human CD11a, also known as LFA-1, and a member of the integrin family of cell adhesion molecules. The antigen is expressed on thymocytes, T and B lymphocytes, and on non-lymphoid cells in the bone marrow. CD11a appears to be an essential component for T-cell activation by antigen-presenting cells and is important for cell-mediated killing. Costimulatory signaling by LFA-1 is achieved by pairing with its ligand, CD54/ICAM-1.

Applications

FC – Quality tested⁵⁻⁸
 IHC-FS – Reported in literature²
 ICC – Reported in literature³
 IP – Reported in literature^{1,4}
 Block – Reported in literature¹

Working Dilutions

Flow Cytometry	FITC and BIOT conjugates	≤ 1 μg/10 ⁶ cells
	PE and PE/CY7 conjugates	≤ 0.2 μg/10 ⁶ cells
For flow cytometry, the suggested use of these reagents is in a final volume of 100 μL		

Other Applications Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.

For Research Use Only. Not for Diagnostic or Therapeutic Use.

Handling and Storage

- The purified (UNLB) antibody is supplied as 0.5 mg of purified immunoglobulin in 1.0 mL of borate buffered saline, pH 8.2. *No preservatives or amine-containing buffer salts added.* Store at 2-8°C.
- The fluorescein (FITC) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The biotin (BIOT) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The R-phycoerythrin (PE) conjugate is supplied as 0.1 mg in 1.0 mL or 0.2 mg in 2.0 mL of PBS/NaN₃ and a stabilizing agent. Store at 2-8°C. **Do not freeze!**
- The low endotoxin, azide-free (LE/AF) antibody is supplied as 0.5 mg of purified immunoglobulin in 1.0 mL of PBS. **Aliquot and store at or below -20°C.**
- The R-phycoerythrin-Cyanine 7 (PE/CY7) conjugate is supplied as 0.1 mg in 1.0 mL of PBS/NaN₃ and a stabilizing agent. Store at 2-8°C. **Do not freeze!**
- Protect fluorochrome-conjugated forms from light. Reagents are stable for the period shown on the label if stored as directed.

Warning

Some reagents contain sodium azide. Please refer to product specific (M)SDS.

References

1. Trowbridge IS, Omary MB. Molecular complexity of leukocyte surface glycoproteins related to the macrophage differentiation antigen Mac-1. *J Exp Med.* 1981;154:1517-24. (Immunogen, IP, Block)
2. Pauschinger M, Noutsias M, Lassner D, Schultheiss HP, Kuehl U. Inflammation, ECG changes and pericardial effusion: whom to biopsy in suspected myocarditis?. *Clin Res Cardiol.* 2006;95:569-83. (IHC-FS, Human Reactivity)
3. Thauland TJ, Koguchi Y, Wetzel SA, Dustin ML, Parker DC. Th1 and Th2 cells form morphologically distinct immunological synapses. *J Immunol.* 2008;181:393-9. (ICC)
4. Lefrançois L, Bevan MJ. Novel antigenic determinants of the T200 glycoprotein expressed preferentially by activated cytotoxic T lymphocytes. *J Immunol.* 1984;135:374-83. (IP)
5. Kato N, Yuzawa Y, Kosugi T, Hobo A, Sato W, Miwa Y, et al. The E-selectin ligand basigin/CD147 is responsible for neutrophil recruitment in renal ischemia/reperfusion. *J Am Soc Nephrol.* 2009;20:1565-76. (FC)
6. Takeshi T, Keisuke N, Takaaki I, Makoto Y, Tatsuji N. Involvement of adhesion molecule in in vitro plaque-like formation of macrophages stimulated with *Aggregatibacter actinomycetemcomitans* lipopolysaccharide. *J Periodont Res.* 2010;45:550-6. (FC)
7. Nagaoka K, Takahara K, Minamino K, Takeda T, Yoshida Y, Inaba K. Expression of C-type lectin, SIGNR3, on subsets of dendritic cells, macrophages, and monocytes. *J Leukoc Biol.* 2010;88:913-24. (FC)
8. Domingos-Pereira S, Decrausaz L, Derré L, Bobst M, Romero P, Schiller JT, et al. Intravaginal TLR agonists increase local vaccine-specific CD8 T cells and human papillomavirus-associated genital-tumor regression in mice. *Mucosal Immunol.* 2013;6:393-404. (FC)

Cy[®] is a registered trademark of GE Healthcare.

TB1555
09-Jan-17

Corporate Offices: 160 Oxmoor Blvd • Birmingham, AL 35209 • USA Mailing Address: P.O. Box 26221 • Birmingham, AL 35260 • USA

Tel: 205.945.1774 • U.S. and Canada: 800.722.2255 • Fax: 205.945.8768

Email: info@southernbiotech.com • Website: www.southernbiotech.com