



## Fc Receptor Saturation Reagent

Cat. No.	Format	Size
5200-01	Purified (UNLB)	50 tests

### Description

<b>Isotype</b>	Rabbit IgG
<b>Source</b>	Normal rabbit serum

### Applications

Quality tested applications include –  
FC

### Working Dilutions

<b>Flow Cytometry</b>	Purified (UNLB) antibody For flow cytometry, the suggested use of these reagents is in a final volume of 100 $\mu$ L	40 $\mu$ L/ $10^6$ cells
<b>Other Applications</b>	Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.	

### Handling and Storage

- The purified (UNLB) IgG is supplied as 50 tests of purified immunoglobulin in 2.0 mL of PBS/ $\text{NaN}_3$ . Store at 2-8°C.
- Reagent is stable for the period shown on the label if stored as directed.

### Warning

Reagent contains sodium azide. Please refer to product specific (M)SDS.

***For Research Use Only. Not for Diagnostic or Therapeutic Use.***

TB5200  
14-Jun-16