


Goat Anti-Mouse IgG_{2b}, Human/Bovine/Horse SP ads

Cat. No.	Format	Size
1093-01	Purified (UNLB)	1.0 mg
1093-02	Fluorescein (FITC)	1.0 mg
1093-04	Alkaline Phosphatase (AP)	1.0 mL
1093-05	Horseradish Peroxidase (HRP)	1.0 mL
1093-08	Biotin (BIOT)	1.0 mg



ELISA plate was coated with purified mouse IgG₁, IgG_{2a}, IgG_{2b}, IgG₃, IgM, and IgA. Immunoglobulins were detected with serially diluted Goat Anti-Mouse IgG_{2b}, Human/Bovine/Horse SP ads-HRP (SB Cat. No. 1093-05).

Description

Specificity	Reacts with the heavy chain of mouse IgG _{2b}
Source	Pooled antisera from goats hyperimmunized with mouse IgG _{2b}
Cross Adsorption	Mouse IgG ₁ , IgG _{2a} , IgG ₃ , IgM, and IgA; human, bovine, and horse serum proteins (SP); may react with immunoglobulins from other species
Purification	Affinity chromatography on mouse IgG _{2b} covalently linked to agarose

Applications

Quality tested applications include –

ELISA
FLISA

Other applications include –

ICC¹
WB¹

Working Dilutions

ELISA	Purified (UNLB) antibody	≤ 1 µg/mL
	AP conjugate	1:2,000 – 1:4,000
	HRP conjugate	1:4,000 – 1:8,000
	BIOT conjugate	1:5,000 – 1:20,000
FLISA	FITC conjugate	1:200 – 1:400
Other Applications	Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.	

For Research Use Only. Not for Diagnostic or Therapeutic Use.

Corporate Offices: 160 Oxmoor Blvd • Birmingham, AL 35209 • USA Mailing Address: P.O. Box 26221 • Birmingham, AL 35260 • USA

Tel: 205.945.1774 • U.S. and Canada: 800.722.2255 • Fax: 205.945.8768

Email: info@southernbiotech.com • Website: www.southernbiotech.com

Handling and Storage

- The purified (UNLB) antibody is supplied as 1.0 mg purified immunoglobulin in 1.0 mL of borate buffered saline, pH 8.2. *No preservatives or amine-containing buffer salts added.* Store at 2-8°C.
- The fluorescein (FITC) conjugate is supplied as 1.0 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The alkaline phosphatase (AP) conjugate is supplied as 1.0 mL in a stock solution of 50 mM Tris/1 mM MgCl₂/50% glycerol, pH 8.0, containing NaN₃ as preservative. Store at 2-8°C or long-term at -20°C.
- The horseradish peroxidase (HRP) conjugate is supplied as 1.0 mL in a stock solution of 50% glycerol/50% PBS, pH 7.4. No preservative added. Store at 2-8°C or long-term at -20°C.
- The biotin (BIOT) conjugate is supplied as 1.0 mg in 2.0 mL of PBS/NaN₃. Store at 2-8°C.
- Protect fluorochrome-conjugated forms from light. Reagents are stable for the period shown on the label if stored as directed.

Warning

Some reagents contain sodium azide. Please refer to product specific SDS.

References

1. SouthernBiotech published data (ICC, WB)