




Hamster Anti-Mouse CD28

Cat. No.	Format	Size
1615-01	Purified (UNLB)	0.5 mg
1615-02	Fluorescein (FITC)	0.5 mg
1615-08	Biotin (BIOT)	0.5 mg
1615-14	Low Endotoxin, Azide-Free (LE/AF)	0.5 mg



BALB/c mouse thymocytes were stained with Hamster Anti-Mouse CD28-FITC (SB Cat. No. 1615-02).

Overview

Clone	PV-1
Isotype	Hamster (Armenian) IgG ₂
Immunogen	C57BL/6N mouse T cell lymphoma EL4 cell line
Specificity	Mouse CD28; Mr 40 kDa
Alternate Name(s)	Tp44, T44

Description

CD28 is a type I disulfide-linked homodimer that is constitutively expressed on most thymocytes, at low density on nearly all CD4⁺ and CD8⁺ peripheral T lymphocytes and at very low levels on NK cells. Its expression is upregulated upon T-cell activation. CD28 is a ligand for CD80 (B7-1) and CD86 (B7-2) on B cells and other antigen presenting cells and plays an important role in the interaction between T cells and B cells. CD28 is a costimulatory receptor involved in many, but not all, T-cell independent immune responses.

Applications

FC – Quality tested ^{1,6}
 IP – Reported in literature ^{1,2}
 Costim – Reported in literature ^{1,3,4}
 Stim – Reported in literature ⁵

Working Dilutions

Flow Cytometry FITC and BIOT conjugates ≤ 3 µg/10⁶ cells
 For flow cytometry, the suggested use of these reagents is in a final volume of 100 µL

Other Applications Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.

For Research Use Only. Not for Diagnostic or Therapeutic Use.

Handling and Storage

- The purified (UNLB) antibody is supplied as 0.5 mg of purified immunoglobulin in 1.0 mL of borate buffered saline, pH 8.2. *No preservatives or amine-containing buffer salts added.* Store at 2-8°C.
- The fluorescein (FITC) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The biotin (BIOT) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The low endotoxin, azide-free (LE/AF) antibody is supplied as 0.5 mg purified immunoglobulin in 1.0 mL of PBS. Contains no preservative; handle under aseptic conditions. Store at 2-8°C or aliquot into smaller volumes and store at -20°C. Avoid multiple freeze / thaw cycles.
- Protect fluorochrome-conjugated forms from light. Reagents are stable for the period shown on the label if stored as directed.

Warning

Some reagents contain sodium azide. Please refer to product specific SDS.

References

1. Abe R, Vandenberghe P, Craighead N, Smoot DS, Lee KP, June CH. Distinct signal transduction in mouse CD4⁺ and CD8⁺ splenic T cells after CD28 receptor ligation. *J Immuno.* 1995;154:985-97. (Immunogen, FC, IP, Costim)
2. Dodson LF, Boomer JS, Deppong CM, Shah DD, Sim J, Bricker TL, et al. Targeted knock-in mice expressing mutations of CD28 reveal an essential pathway for costimulation. *Mol Cell Biol.* 2009;29:3710-21. (IP)
3. Kanzaki M, Wada J, Sugiyama K, Nakatsuka A, Teshigawara S, Murakami K, et al. Galectin-9 and T cell immunoglobulin mucin-3 pathway is a therapeutic target for type 1 diabetes. *Endocrinology.* 2012;153:612-20. (Costim)
4. Tomimori Y, Tanaka Y, Goto M, Fukuda Y. Repeated topical challenge with chemical antigen elicits sustained dermatitis in NC/Nga mice in specific-pathogen-free condition. *J Invest Dermatol.* 2005;124:119-24. (Costim)
5. Zhang T, Fresnay S, Welty E, Sangrampurkar N, Rybak E, Zhou H, et al. Selective CD28 blockade attenuates acute and chronic rejection of murine cardiac allografts in a CTLA-4-dependent manner. *Am J Transplant.* 2011;11:1599-1609. (Stim)
6. Jacobsen J, Haabeth OW, Tveita AA, Schjetne KW, Munthe LA, Bogen B. Naive idiotope-specific B and T cells collaborate efficiently in the absence of dendritic cells. *J Immunol.* 2014;192:4174-83. (FC)

Licensed from US Navy under US Patent No. 5,948,893.

TB1615
08-Oct-21