SouthernBiotech



Goat Anti-Rat IgG

Cat. No.	Format	Size
3030-01	Purified (UNLB)	1.0 mg
3030-02	Fluorescein (FITC)	1.0 mg
3030-03	Rhodamine (TRITC)	1.0 mg
3030-04	Alkaline Phosphatase (AP)	1.0 mL
3030-05	Horseradish Peroxidase (HRP)	1.0 mL
3030-06	β-galactosidase (BGAL)	1.0 mL
3030-08	Biotin (BIOT)	1.0 mg
3030-09	R-phycoerythrin (PE)	0.5 mg
3030-30	Alexa Fluor [®] 488 (AF488)	1.0 mg
3030-31	Alexa Fluor [®] 647 (AF647)	1.0 mg
3030-32	Alexa Fluor [®] 555 (AF555)	1.0 mg



FLISA plate was coated with purified rat IgG and IgM. Immunoglobulins were detected with serially diluted Goat Anti-Rat IgG-PE (SB Cat. No. 3030-09).

Description

Specificity	Reacts with the heavy chain of rat IgG
Source	Pooled antisera from goats hyperimmunized with rat IgG
Cross Adsorption	Rat IgM; may react with immunoglobulins from other species
Purification	Affinity chromatography on pooled rat IgG covalently linked to agarose

Applications

Quality tested applications include – ELISA $^{1\mbox{-}5}$

FLISA

Other referenced applications include -

FC ^{6,7} IHC-FS ^{8,9} ICC ^{10,11} WB ¹²⁻¹⁵

Working Dilutions

ELISA	AP conjugate HRP conjugate BGAL conjugate BIOT conjugate	1:2,000 - 1:4,000 1:4,000 - 1:8,000 1:500 1:5,000 - 1:20,000
FLISA	FITC, TRITC, AF488, and AF555 conjugates PE and AF647 conjugates	1:100 – 1:400 ≤ 1 μg/mL
Other Applications	Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.	

For Research Use Only. Not for Diagnostic or Therapeutic Use.

Handling and Storage

- The purified (UNLB) antibody is supplied as 1.0 mg purified immunoglobulin in 1.0 mL of borate buffered saline, pH 8.2. *No preservatives or amine-containing buffer salts added*. Store at 2-8°C.
- The fluorescein (FITC), rhodamine (TRITC), Alexa Fluor[®] 488 (AF488), Alexa Fluor[®] 555 (AF555), and Alexa Fluor[®] 647 (AF647) conjugates are supplied as 1.0 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The alkaline phosphatase (AP) conjugate is supplied as 1.0 mL in a stock solution of 50 mM Tris/1 mM MgCl₂/50% glycerol, pH 8.0, containing NaN₃ as preservative. Store at 2-8°C or long-term at -20°C.
- The horseradish peroxidase (HRP) conjugate is supplied as 1.0 mL in a stock solution of 50% glycerol/50% PBS, pH 7.4. No preservative added. Store at 2-8°C or long-term at -20°C.
- The β-galactosidase (BGAL) conjugate is supplied as 1.0 mL of stock solution in 50% glycerol/50% PBS containing NaN₃ as preservative. Store at 2-8°C or long-term at -20°C.
- The biotin (BIOT) conjugate is supplied as 1.0 mg in 2.0 mL of PBS/NaN₃. Store at 2-8°C.
- The R-phycoerythrin (PE) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN₃ and a stabilizing agent. Store at 2-8°C. **Do not** freeze!
- Protect fluorochrome-conjugated forms from light. Reagents are stable for the period shown on the label if stored as directed.

Warning

Some reagents contain sodium azide. Please refer to product specific SDS.

References

- 1. Laumonier T, Mohacsi PJ, Matozan KM, Banz Y, Haeberli A, Korchagina EY, et al. Endothelial cell protection by dextran sulfate: a novel strategy to prevent acute vascular rejection in xenotransplantation. Am J Transplant. 2003;4:181-7. (ELISA)
- 2. Wernhoff P, Olofsson P, Holmdahl R. The genetic control of rheumatoid factor production in a rat model of rheumatoid arthritis. Arthritis Rheum. 2003;48:3584-96. (ELISA)
- 3. Eto Y, Yoshioka Y, Ishida T, Yao X, Morishige T, Narimatsu S, et al. Optimized PEGylated adenovirus vector reduces the anti-vector humoral immune response against adenovirus and induces a therapeutic effect against metastatic lung cancer. Biol Pharm Bull. 2010;33:1540-4. (ELISA)
- 4. Cheong C, Matos I, Choi J, Schauer JD, Dandamudi DB, Shrestha E, et al. New monoclonal anti-mouse DC-SIGN antibodies reactive with acetonefixed cells. J Immunol Methods. 2010;360:66-75. (ELISA)
- 5. Zuest R, Valdes I, Skibinski D, Lin Y, Toh YX, Chan K, et al. Tetravalent dengue DIIIC protein together with alum and ODN elicits a Th1 response and neutralizing antibodies in mice. Vaccine. 2015;33:1474-82. (ELISA)
- 6. Achour TN, Sentis S, Teyssier C, Philippat A, Lucas A, Corbo L, et al. Transcriptional repression of estrogen receptor α signaling by SENP2 in breast cancer cells. Mol Endocrinol. 2014;28:183-96. (FC)
- 7. Scheibe M, Arnoult N, Kappei D, Buchholz F, Decottignies A, Butter F, et al. Quantitative interaction screen of telomeric repeat-containing RNA reveals novel TERRA regulators. Genome Res. 2013;23:2149-57. (FC)
- 8. Huizinga R, Easton AS, Donachie AM, Guthrie J, van Rijs W, Heikema A, et al. Sialylation of Campylobacter jejuni lipo-oligosaccharides: impact on phagocytosis and cytokine production in mice. PLoS One. 2012;7(3):e34416. (IHC-FS)
- 9. Duehrkop C, Banz Y, Spirig R, Miescher S, Nolte MW, Spycher M, et al. C1 esterase inhibitor reduces lower extremity ischemia/reperfusion injury and associated lung damage. PLoS One. 2013;8(8):e72059. (IHC-FS)
- 10. Alliot F, Lecain E, Grima B, Pessac B. Microglial progenitors with a high proliferative potential in the embryonic and adult mouse brain. Proc Natl Acad Sci USA. 1991;88:1541-5. (ICC)
- 11. Taoudi S, Morrison AM, Inoue H, Gribi R, Ure J, Medvinsky A. Progressive divergence of definitive haematopoietic stem cells from the endothelial compartment does not depend on contact with the foetal liver. Development. 2005;132:4179-91. (ICC)
- 12. O'Reilly LA, Print C, Hausmann G, Moriishi K, Cory S, Huang DC, et al. Tissue expression and subcellular localization of the pro-survival molecule Bcl-w. Cell Death Differ. 2001;8:486-94. (WB)
- 13. O'Reilly LA, Ekert P, Harvey N, Marsden V, Cullen L, Vaux DL, et al. Caspase-2 is not required for thymocyte or neuronal apoptosis even though cleavage of caspase-2 is dependent on both Apaf-1 and caspase-9. Cell Death Differ. 2002;9:832-41. (WB)
- 14. Armitage CW, O'Meara CP, Harvie MC, Timms P, Blumberg RS, Beagley KW. Divergent outcomes following transcytosis of IgG targeting intracellular and extracellular chlamydial antigens. Immunol Cell Biol. 2014;92:417-26. (WB)
- 15. Lou Y, Wu C, Wu M, Xie C, Ren L. The changes of neutrophil gelatinase-associated lipocalin in plasma and its expression in adipose tissue in pregnant women with gestational diabetes. Diabetes Res Clin Pract. 2014;104:136-42. (WB)

Alexa Fluor[®] 488, 647, and 555 are provided under an Intellectual property license from Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For information on purchasing a license to this product for any other use, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com.